





DP Hoof Photo

10:54:30

Owner: SampleOwner Animal: SampleHorse Date: 3-Nov-2022



The 'Splay' measurement helps to assess the frontal view of the hoof. Horses with low splay values tend to be contracted or clubby, and horses with high splay tend to be weak walled and/or have weak bars. No one value is 'correct' for all animals, but extreme values are to be avoided.



The Splay for this hoof is 121.4 percent. This value is indicated on the histogram to the left by the vertical black line labelled 'This Hoof'. Splay values substantially higher or lower than normal are to be avoided, if possible.

The Splay value is in the 1st (low) quartile when compared to a large group horses which have been measured.

Data from 1,880 hooves of mixed breeds.